

Implant compatibility chart



Manufacturer	Implant	Titanium	GoldHue	Zirconia, white	Zirconia, shaded	Screwdriver type
Astra Tech	Internal					
	OsseoSpeed™/OsseoSpeed™ TX 3.0 - Yellow	x	x			0.050" Hex screwdriver
	OsseoSpeed™/OsseoSpeed™ TX 3.5, 4.0 - Aqua	x	x	x	x	
	OsseoSpeed™/OsseoSpeed™ TX 4.5, 5.0 - Lilac	x	x	x	x	
OsseoSpeed™ 4.0 - Blue (Previous)	x	x				
BioHorizons	Internal					
	Internal 3.5, 4.5, 5.7	x	x	x	x	1.25 mm Hex screwdriver
	Tapered Internal 3.5, 4.5, 5.7	x	x	x	x	
	External					
	External Implant 3.5	x	x			0.050" Hex screwdriver
	External Implant 4.0, 5.0	x	x	x	x	
Biomet 3i	Internal					
	Certain™ 3.25, 4.0, 5.0, 6.0	x	x	x	x	0.048" Hex screwdriver
	Certain™ 4/3, 5/4	x	x	x	x	
	Certain™ Prevail™ 3/4/3, 4/5/4, 5/6/5	x	x	x	x	
	Certain™ Prevail™ 4/3, 5/4, 6/5	x	x	x	x	
	Certain™ XP 4/5, 5/6	x	x	x	x	
	External					
	MicroMiniplant™ 3.25	x	x			0.048" Hex screwdriver
	Miniplant™ 3.25	x	x	x	x	
	Osseofite™ 3.75, 4.0, 5.0, 6.0	x	x	x	x	
Osseofite™ XP 4/5, 5/6	x	x	x	x		
Dentsply Friadent	Internal					
	XiVE S 3.0	x	x			0.048" Hex screwdriver
XiVE S/Frialit 3.4, 3.8, 4.5, 5.5	x	x	x	x		
Keystone Dental	External					
	Restore™ Small Diameter 3.3, 3.4	x	x			0.048" Hex screwdriver
	Restore™ Regular Diameter 3.75, 4.0	x	x	x	x	
Restore™ Wide Diameter 5.0, 5.5, 6.0	x	x				
Nobel Biocare	Internal					
	NobelActive™ NP 3.5	x	x	x	x	Nobel Biocare Unigrip
	NobelActive™ RP 4.3, 5.0	x	x	x	x	
	NobelReplace™ NP 3.5	x	x	x	x	
	NobelReplace™ RP 4.0	x	x	x	x	
	NobelReplace™ WP 5.0	x	x	x	x	
	NobelReplace™ 6.0	x	x	x	x	
	External					
	Brånemark System™ NP 3.3	x	x			Nobel Biocare Unigrip
	Brånemark System™ RP 3.75, 4.0	x	x	x	x	
	Brånemark System™ WP 5.0, 5.5	x	x	x	x	
	SteriOss Replace™ 3.1, 3.5, 4.3, 5.0, 6.0	x	x			0.050" Hex screwdriver
SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0	x	x				
Straumann	Internal					
	Bone Level 3.3 NC (Narrow CrossFit™)	x	x			Straumann screwdriver compatible
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	x	x	x	x	
	Standard/Standard Plus 4.8 RN (Regular Neck)	x	x			
External						
Standard Plus 3.3 NN (Narrow Neck)	x	x			Straumann screwdriver compatible	
Sybron Implant Solutions	External					
	Endopore 3.5, 4.1, 5.0 - external	x	x			0.050" Hex screwdriver
Entegra 3.5, 4.1, 5.0 - external	x	x				
Zimmer Dental	Internal					
	Tapered Screw-Vent™ 3.5, 4.5, 5.7	x	x	x	x	1.25 mm Hex screwdriver
	Screw-Vent™ 3.3, 3.7, 4.5	x	x	x	x	
SwissPlus™ 4.8	x	x				

Turnaround times

Turnaround times (in-house business days) for all case sizes		
Titanium	Gold-shaded titanium	Zirconia
2	4	4

Note: Turnaround time listed begins the first business day after receipt of all case information and material required and does not include shipping, which is determined by your choice of delivery method. Turnaround for lab-based scanned cases is one day less than the times shown.

Turnaround times are based on the following:

- Design approvals to process case must be granted by 10 am your local time.
- Lab-scanned cases, including rescans, must be uploaded to Atlantis™ WebOrder by 3 pm (EST) in order to save a day of processing compared with the times listed above.
- Material, on-hold and case information issues must be resolved by 8 pm (EST).

Other services

Service	Additional turnaround days
Laboratory services	add 1 day
Provisional restoration	add 5 days
Removable diagnostic preparation	add 2 days

- For the most current implant compatibility, please visit www.atlantisabutment.com or contact Customer Service at 800-531-3481.
- All orders for Atlantis™ Abutments are shipped with the corresponding abutment screw needed at no extra charge.

Note: The implant manufacturer's recommended torque should always be used when placing the abutment screw.

